

# BWT bestmin


## Special filter system for water remineralisation downstream of reverse osmosis systems



BWT bestmin COFFEE is an innovative filter that is ideally matched to **BWT bestaqua reverse osmosis systems**. Reverse osmosis allows dissolved substances capable of leading to undesirable deposits – such as limescale and gypsum – to be removed from the water. The process thereby also eliminates chloride, protecting coffee machines from corrosion. Reverse osmosis thus produces an exceptionally pure water, albeit water that lacks something essential: minerals! They need to be added using the BWT bestmin COFFEE filter in order to perfect the brewing process and the sensory quality of the coffee.


In this filter, special mineral granules ensure **targeted remineralisation**, while colloids and bacteria (log6) are trapped and removed from the water in a membrane module (0,1 µm), enhancing the purity and microbiological integrity of the water.

Reverse osmosis systems in the BWT bestaqua range and BWT bestmin COFFEE for remineralisation are an **ideal combination**, marrying protection of the machine with the highest quality of coffee, regardless of the local water quality and even if using water with a high salt content. BWT bestmin COFFEE filters are food-safe and can be used downstream of any quality, and correctly configured reverse osmosis system and with any quality of drinking water worldwide.



## Technical data


BWT bestmin COFFEE 100	
Typical capacity in litres	10,000
Height excluding handle (A) in mm	475
Height with handle (B) in mm	500
Connection height (C) in mm	425
Ø filter candle (D) in mm	130
Weight in kg (dry/wet)	4.0/4.5
Filter candle order no.	FS24M99A00



### Operating conditions

Feed pressure:	2–8 bar
Water temperature:	4–30 °C
Ambient temperature:	4–40 °C

## BWT bestaqua COFFEE – The complete system



## BWT bestmin COFFEE – The mineral-giver

- For perfect remineralisation of water following reverse osmosis treatment
- Setting the mineral balance of Mg and Ca to an ideal ratio
- Neutral pH value for optimum coffee extraction
- With additional membrane filter stage for removing suspended solids (colloids) and microorganisms



## Accessories

Description	Order number
<b>BWT besthead FLEX with 2 90° elbow connectors, FLEX Insert, M 3/8"</b>	FS00Z20A00
<b>BWT besthead FLEX 1 connection set</b>	
– BWT besthead FLEX with 2 90° elbow connectors, FLEX Insert, M 3/8"	
– Connector hose DN 8, 1.5 m, with ball valve, FLEX Insert 90° elbow, cap nut FM 3/8"	
– Straight connector, FLEX Insert, M 3/8"	FS00Z39A00
– 90° elbow connector, FLEX Insert, M 3/8"	
– Adapter FM 3/4" x M 3/8"	
<b>BWT Test Kit for hot beverages</b>	YS00Y01A00



Subject to change without notice.

Manufacturer:  
**BWT water + more GmbH**  
 Walter-Simmer-Str. 4 · A-5310 Mondsee Austria  
 Tel. +43 6232 5011-1164

Sales:  
**BWT water+more Deutschland:** info@water-and-more.de · Tel.: +49 611 58019-0  
**BWT water + more Austria:** kundenservice.wm@bwt-group.com · Tel.: +43 6232 5011-1164  
**BWT Aqua AG Schweiz:** info@bwt-aqua.ch · Tel.: +41 61 755 88 99

[www.bwt-wam.com](http://www.bwt-wam.com)

**BWT**  
 water + more